

10011742 The Alaris® Tri-Site Thermometer model 10011742 is a portable single probe instrument equipped with a pulse timer. It offers both predictive and continuous oral and axillary measurements. All Alaris® Tri-Site Thermometers feature a three-button front panel for easy access to functions. Comes with three AA batteries an anti-theft security base, carrying strap, backlight, probe cover dispenser, and a three-year limited warranty.

10011743 The Alaris® Tri-Site Thermometer model 10011743 is a portable single probe instrument equipped with a pulse timer. It offers both predictive and continuous rectal measurements. All Alaris® Tri-Site Thermometers feature a three-button front panel for easy access to functions. Comes with three AA batteries, anti-theft security base, carrying strap, backlight, probe cover dispenser, and a three-year limited warranty.

10011905 The Alaris® Tri-Site Thermometer model 10011905 is a wall mount dual probe instrument equipped with a pulse timer. It offers both predictive and continuous oral, rectal and axillary measurements. All Alaris® Tri-Site Thermometers feature a three-button front panel for easy access to functions. Comes with a security base designed for attachment to a wall or other secure location. The 1011905 comes with three AA batteries, an anti-theft security base, backlight, probe cover dispenser, and a three-year limited warranty.


10011906 The Alaris® Tri-Site Thermometer model 10011906 is a portable single probe instrument equipped with a pulse timer. It offers both predictive and continuous oral and axillary measurements. All Alaris® Tri-Site Thermometers feature a three-button front panel for easy access to functions. Model 1001906 comes with security base designed for attachment to a wall or other secure location. Comes with three AA batteries, an anti-theft security base, backlight, probe cover dispenser, and three-year limited warranty.

10011907 The Alaris® Tri-Site Thermometer model 10011907 is a portable single probe instrument equipped with a pulse timer. It offers both predictive and continuous rectal measurements. All Alaris® Tri-Site Thermometers feature a three-button front panel for easy access to functions. Model 10011907 comes with a security base designed for attachment to a wall or other secure location. The 10011907 come with three AA batteries, an anti-theft security base, backlight, probe cover dispenser, and a three-year limited warranty.


2031A Thermometer stand for Turbo Temp® Thermometers and Alaris® Tri-Site Thermometers.



2117A Home base unit for use with portable Alaris® Tri-Site and Turbo Temp® Thermometers.



2118A Home base unit for use with wall mount Alaris® Tri-Site and Turbo Temp® Thermometers.



3887 Blue oral and probe with long white cord (80") for use with Alaris® Tri-Site Thermometer models 1001741, 1001742, 1001905, and 1001906. (Uses P-850A probe covers.)

4888 Red rectal probe with long white cord (80") for use with Alaris® Tri-Site Thermometer models 1001741, 1001743, 1001905, and 1001907. (Uses P-850A probe covers.)

P-850A For use with Temp Plus® II, Turbo Temp®, Alaris® Tri-Site Thermometer and vital Check® vital signs monitor series 4400/4500. GE Dinamap monitors with Temp Plus® II of Turbo Temp® Thermometers also use these probe covers. Count: 20/box, 10 boxes/dispensers, 25 dispensers/case.



TE1811 Probe simulator for Turbo Temp® and Alaris® Tri-Site Thermometers.