



- Alaris® Products home
- **Products**
- + Infusion systems
- Guardrails® software
- + Networking products
- + Safety disposables
- + Thermometry
- + Services and solutions
- + OEM
- + Support




Home > Products and services > Alaris® Products > Products > e-Catalog >


e-Catalog search results

-->


2031A Thermometer stand for Turbo Temp® Thermometers and Alaris® Tri-Site Thermometers.



2117A Home base unit for use with portable Alaris® Tri-Site and Turbo Temp® Thermometers.




2118A Home base unit for use with wall mount Alaris® Tri-Site and Turbo Temp® Thermometers.




2180C The Turbo Temp® Thermometer model 2180C is a portable single probe instrument equipped with a pulse timer. It offers predictive oral capability, as well as oral/ axillary site measurements in the monitoring setting. Includes three AA batteries an anti-theft security base, carrying strap, backlight, probe cover dispenser, and a three-year limited warranty.







2180CX01EE The Turbo Temp® Thermometer model 2180CX01EE is a portable dual probe instrument equipped with a pulse timer. It offers predictive oral and rectal capabilities as well as oral, rectal and axillary measurements in the monitoring mode. The 2180C comes with three AA batteries, an anti-theft security base, carrying strap, backlight, probe cover dispenser, and a three-year limited warranty.

2185B The Turbo Temp® Thermometer model 2185B is a wall mount, single probe instrument equipped with a pulse timer. It offers predictive oral capability as well as oral and axillary measurements in monitoring mode. Model 2185B comes with a security base designed for attachment to a wall or other secure location. It comes with three AA batteries, an anti-theft security base, backlight, probe cover dispenser, and a three-year limited warranty.







2185BX01EE The Turbo Temp® Thermometer model 2185BX01EE is a wall mount dual probe instrument equipped with a pulse timer. It offers predictive oral and rectal capabilities as well as oral, rectal and axillary measurements in the monitor mode. Model 2185BX01EE comes with security base design for attachment to a wall or other secure location. The 2185BX01EE comes with three AA batteries, an anti-theft security base, backlight, probe cover dispenser, and a three-year limited warranty.




2885C Blue oral probe with short gray cord (40") for Turbo Temp® Thermometer Model 2180C (uses P850A probe covers).



2886C Red rectal probe with short gray cord (40") for Turbo Temp® Thermometer Model 2180. (Uses P-850A Probe Covers).



2887B Blue oral probe with long gray cord (80") for use with Turbo Temp® Thermometer Model 2185 (uses P-850A probe covers).



2888B Red rectal probe with long gray cord (80") for use with Turbo Temp® Thermometer Model 2185 (uses P-850A probe covers.)

P-850A For use with Temp Plus® II, Turbo Temp®, Alaris® Tri-Site Thermometer and vital Check® vital signs monitor series 4400/4500. GE Dinamap monitors with Temp Plus® II of Turbo Temp® Thermometers also use these probe covers. Count: 20/box, 10 boxes/dispensers, 25 dispensers/case.



TE1811 Probe simulator for Turbo Temp® and Alaris® Tri-Site Thermometers.

[back to top](#)

<!-- ----- * END: DISPLAY PRODUCTS
----- --> **Please note: These products are only available in the USA and Canada.**

Feedback: [Alaris® Products](#)

[Home](#) | [About us](#) | [Careers](#) | [Investors](#) | [Media](#) | [Contact us](#) | [Legal](#) | [Privacy policy](#)