## BCI™ 3180 PULSE OXIMETER

## SPECIFICATIONS (Nominal):

Displays:

Display Backlit LCD
Resolution 240 X 64 pixels
Waveform Plethysmogram
Graphic Trends SpO<sub>2</sub>, Pulse Rate,

Pulse Signal Strength

 $SpO_2$ 

Pulse Rate

Pulse Strength Logarithmically scaled

8-segment bar graph

Indicators:

Low Battery Low battery icon when 10 minutes of

battery life remain

Audio Status Yellow LED and displayed symbol

indicates alarm audio paused

(2 min.) or off

SpO<sub>2</sub> Sensor Message indicates sensor problem

SpO<sub>2</sub> Search Message indicates monitor is searching

for a pulse

Small Pulse Message indicates oximeter signal

strength is low

Charge Green CHARGE LED indicates external

power and battery is charging

Lost Pulse Message indicates monitor has

lost the signal

Sp0<sub>2</sub>:

Range 0-100%

Accuracy

Adult ±2% @70-100% SpO<sub>2</sub>
Neonate ±3% @70-100% SpO<sub>2</sub>
<70% unspecified

Alarm Range High=50-100% & OFF (1% steps)

Low=50-99% & OFF (1% steps)

Averaging 4, 8, or 16 second average

## SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE

Please read the Instructions for Use supplied with the product for detailed instructions, warnings, and cautions.

Pulse Rate:

Range 30-254 bpm

Accuracy ±2 bpm or ±2% (whichever is greater)

@ 30-254 bpm

Alarm Range High=30-254 bpm & OFF

Low=30-254 bpm & OFF

Averaging 8 or 16 second average

Dimensions:

8"W x 6.5"H x 5"D

(20.3cm x 16.5cm x 12.7cm)

Weight:

4.5lb (2kg)

AC Power Requirements:

Two voltage ranges selected by external switch

90-125VAC 50/60Hz 200-240VAC 50/60Hz

Battery Life:

Standard lead-acid 6 VDC battery:

Internal rechargeable
4 hour continuous use (new)

Printer Output:

Serial infrared link-real time patient data

RS232 trend table output

Environmental Specifications:

Operating Temperature 32°F to 104°F (0°C to 55°C)
Storage Temperature -40°F to 167°F (-40° to 75°C)

Relative Humidity 15-95% (operating) (non condensing) 10-95% (storage)

Standards:

EN60601-1 (IEC 601-1)
UL 2601-1 and CAN/CSA 601.1

CAN/CSA NRTL/C

Regulatory Compliance:

510(k): K020350





ISO 9001 / EN 46001

## Catalog Number: 3180

**Q-Line for U.S.A. Customers** call 1-800-348-1842 for emergency and clinical application information 24 hours a day, 365 days a year. Experienced Clinical Specialists will return your call and assist with any questions you may have.

For more information or for the name of your representative, contact BCI, Inc.



BCI, Inc., N7 W22025 Johnson Road, Waukesha, Wisconsin 53186 USA Phone: 262-542-3100 Fax: 262-542-0718 Toll-Free: 800-558-2345 e-mail: info@smiths-bci.com Internet: http://www.smiths-bci.com

smiths

Smiths Medical - a part of Smiths Group plc

