

Fast electronic thermometer system developed for patient safety



- Dedicated probe storage for infection control
- Fast and accurate: oral temperatures in approximately 7 seconds*
- Large, easy-to-read screen
- Ambient light sensor
- Photosensitive backlight
- Pulse timer
- Anti-theft protection

Single Probe
Oral Thermometer



The TurboTemp® Thermometer was developed with special attention to patient safety, this innovative "dedicated probe" electronic thermometer offers many important safety and convenience features.



The easy-to-use choice for safety

Fast and accurate

Oral temperatures in approximately 7 seconds*

Focused on safety

Oral and rectal probe thermometers to reduce the risk of cross contamination

Convenience

Full 60-second pulse timer. Automatic back-lit screen for night time readings

Productivity

Temperature recall at the press of a button

Anti-theft

Electronic loss prevention deters theft and emits audible signal if removed from home base for an extended period

Specifications

Models 2180 and 2185

Normal operation	96°F to 106.0°F 35.6°C to 41.1°C
Monitor mode	80°F to 106.0°F 26.7°C to 41.1°C
Pulse/Respiration timer	Audible signals at 15, 30 and 60 sec.
Visual display	every 5 sec.
Accuracy	±0.2°F (monitor mode when tested in a ±0.1°C calibrated water bath)
Battery	3 AA size alkaline cells

Dimensions

Length	17.5 cm. 6.88 in.
Width	7.30 cm. 2.88 in.
Height	5.7 cm. 2.23 in.
Weight	.32 kg./ .7 lb. (.9 lb. w/batteries)
ALARIS® Probe Cover 3 3/4 in. long	Model P850 Packaged 20 per box, 200 per carton 5,000 per shipping box

*Refer to directions for use



Single Probe Rectal Thermometer



Economical ALARIS® Probe Covers for patient safety and accuracy

For more information, contact your ALARIS Medical Systems Sales Consultant at 1-800-482-4822, in Canada 1-800-387-8309, Fax 1-858-458-7760, or visit our web site at www.alarismed.com

ALARIS®
MEDICAL SYSTEMS
Medication Safety
at the Point of Care™

ALARIS Medical Systems, Inc.
Worldwide Headquarters,
10221 Wateridge Circle,
San Diego, California 92121-2772

Fax: 1-858-458-7760
Customer Service: 1-800-482-4822
Website: www.alarismed.com