

KANGAROO* 224

Enteral Feeding Pump



tyco

Healthcare

KENDALL

KANGAROO* 224

Enteral Feeding Pump

Rotary peristaltic pump with very simple operation and exceptionally low maintenance

The KANGAROO* 224 is a rotary peristaltic pump for the administration of probe feeds. It features exceptional ease of use: Power-on and power-off at the touch of a button on the keypad, and setting of the volume to be administered via the rotary knob. With 3 simple controls, the KANGAROO* 224 is ideal for home-fed patients.



KANGAROO* 224

REF	Case Quantity	Model
8884332401	1	European version (two pin plug)
8884332400	1	English version (three pin plug)

Mobile infusion stand available as an accessory



PRODUCT FEATURES:

- Continuous delivery: 5 - 300 ml/hr
- Closing pressure: 0.8 bar
- Battery life: 24 hrs at average flow rate
- Battery charging time: approx. 15 hrs
- Audible and visual alarm system

K224 Accessories

REF	Description	Case Quantity
1814-340228	Walking stick dripstand with wheels + long hook	1
1814-340226**	Walking stick dripstand with wheels + short hook	1
1814-340210**	Walking stick dripstand without wheels	1
1814-340236**	Table top dripstand	1



Easy and safe:
Unique alarm alert for an improperly loaded set.

** without CE

* KANGAROO is a trademark of Sherwood Services AG. © 2003 Tyco Healthcare Group LP. All rights reserved.

P-EF-K224/GB

tyco
Healthcare
KENDALL

Tyco Healthcare UK Limited
154 Fareham Road
GOSPORT Hampshire
PO 13 0AS - UK
Tel: +44 (0) 1329 224226
Fax: +44 (0) 1329 224334
www.tycohealthcare.co.uk

Tyco Healthcare
Regional Export
154 Fareham Road
GOSPORT Hampshire
PO 13 0AS
UK
Tel: +44 (0) 1329 224258
Fax: +44 (0) 1329 224086

Tyco Healthcare Nederland BV
Hogeweg 105
NL-5301 LL ZALTBOMMEL
The NETHERLANDS
Tel: +31 418 57 66 00
Fax: +31 418 57 67 91
TycoHealthNed@TycoHealth.com
www.tycohealthcare.nl

Tyco Healthcare
Birchwood Court Montrose Road
P.O. Box 8108, Halfway House
1685 MIDRAND
SOUTH AFRICA
Tel: +27 (0) 11 265 2000
Fax: +27 (0) 11 265 2050

Tyco Healthcare Belgium N.V.S.A.
Generaal de Wittelaan 9/5
B-2800 MECHELEN
BELGIUM
Tel: +32 (0) 15 29 44 50
Fax: +32 (0) 15 29 44 55
tycohealthcare.Belgium@tyco-healthcare.com

Tyco Healthcare Norden AB
Gårdsvägen 18, Box 710
SE-169 27 Solna
SWEDEN
Tel: +46-8 585 605 00
Kundservice +46-8 585 605 05
Fax: +46-8 585 605 06

Tyco Healthcare Denmark
Langebrogade 6E, 4.th.
DK-1411 København K
DENMARK
Tel: +45-70 27 52 50
Kundeservice +45-70 27 52 50
Fax: +45-70 27 51 50

Tyco Healthcare Norge AS
Holmengaten 24
N-1394 Nesbru
NORWAY
Tel: +47-66 85 22 00
Fax: +47-66 85 22 11

Tyco Healthcare Finland
Köydenpunojankatu 2D, 6 krs, PL340
FIN-00181 Helsinki
FINLAND
Puh +358-(0)9-6226 840
Asiakaspalvelu
+358-(0)9-6226 8418
Fax: +358-(0)9-6226 8460