


# Tracheostomy



**PORTEX**<sup>TM</sup>

Portex® (Polyvinyl Chloride – PVC)

Bivona® (Silicone)

Tracheostoma Vents

Percutaneous

Tracheostomy Kits

Cricothyrotomy Kits

Tracheostomy Ancillary

Products and Accessories


For More  
Information  
Call Customer  
Service at  
800.258.5361

## Tracheostomy

<b>01 Percutaneous Tracheostomy Kits</b>
p2 ULTRAPerc®
p2 Per-fit™
<b>02 Portex® Tracheostomy Tubes (Polyvinyl Chloride – PVC)</b>
p7 Adult Tracheostomy Tubes
p7 Disposable Inner Cannula (D.I.C.®) Tubes
p8 D.I.C.® Laryngectomy Tubes
p9 Portex® Blue Line ULTRA® Tracheostomy Tubes
p11 Portex® Blue Line® Tracheostomy Tubes
p11 Talk Tube
p12 Adjustable and Extra Length Tracheostomy Tubes
p13 Neonatal and Pediatric Tracheostomy Tubes
<b>03 Bivona® Tracheostomy Tubes (Silicone)</b>
p17 Adult Tracheostomy Tubes
p17 Adjustable and Extra Length Tracheostomy Tubes
p17 Silicone Tracheostomy Tubes
p19 Talk Tubes
p19 Neonatal and Pediatric Tracheostomy Tubes
p19 Uncuffed Neonatal and Pediatric Tracheostomy Tubes
p20 Cuffed Neonatal and Pediatric Tracheostomy Tubes
p21 Specialty Pediatric Tracheostomy Tube
p22 Customized Tracheostomy Tube Services
<b>04 Bivona® Tracheostoma Vents</b>
p25 Tracheostoma Vents
<b>05 Cricothyrotomy Kits</b>
p27 Mini-Trach® II Cricothyrotomy Kit
p27 Nu-Trake® Adult Emergency Cricothyrotomy Kit
<b>06 Tracheostomy Ancillary Products and Accessories</b>
p30 Velcro Tracheostomy Tube Holders
p30 Tracheostomy Tube Straps
p30 Tracheostomy Care Trays
p30 Tracheostomy Educational Aids
p31 Tracheostomy Accessories

# Tracheostomy

## Percutaneous Tracheostomy Kits



<b>ULTRAPerc® (Single-Stage Dilation)</b>	<b>p2</b>
Portex® ULTRAPerc® Percutaneous Dilation Tracheostomy Kits with Blue Line ULTRA® Tube	p2
Portex® ULTRAPerc® Percutaneous Dilation Tracheostomy Kits without Tube	p2
Portex® ULTRAPerc® Local Infiltration Tray	p2

<b>Per-fit™ (Serial Dilation)</b>	<b>p2</b>
Per-fit™ Percutaneous Tracheostomy Kits	p2
Replacement Per-fit™ Percutaneous Tracheostomy Tubes	p2
Inner Cannula for Per-fit™ Percutaneous Tracheostomy Tubes	p2



100/591/080

## Portex® ULTRAPerc® Percutaneous Dilation Tracheostomy Kits with Blue Line ULTRA® Tube

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
100/591/070	7.0 mm	7.0 mm	10.5 mm	70.0 mm	1
100/591/080	8.0 mm	8.0 mm	11.9 mm	75.5 mm	1
100/591/090	9.0 mm	9.0 mm	13.3 mm	81.0 mm	1

Includes two inner cannula, curved introducer with handle, disconnection wedge, tracheostomy tube holder, cleaning brush, prep and procedural components.



100/572/000

## Portex® ULTRAPerc® Percutaneous Dilation Tracheostomy Kits without Tube

PRODUCT CODE	DESCRIPTION	UNITS PER CASE
100/572/000	Includes 21fr, 24fr, 27fr Obturators, Plus Prep and Procedural Components	1



100/572/LIT

## Portex® ULTRAPerc® Local Infiltration Tray

PRODUCT CODE	DESCRIPTION	UNITS PER CASE
100/572/LIT	Compact Double Tray System with Drugs, Prep and Procedural Components	10



535090

## Per-fit™ Percutaneous Tracheostomy Kits

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	COLOR CODE	UNITS PER CASE
535070	7.0 mm	7.0 mm	9.6 mm	82.0 mm	Green	2
535080	8.0 mm	8.0 mm	10.9 mm	86.0 mm	White	2
535090	9.0 mm	9.0 mm	12.3 mm	93.0 mm	Blue	2

Each kit includes four size-specific dilators and obturator / dilator plus prep and procedural components.



536090

## Replacement Per-fit™ Percutaneous Tracheostomy Tubes

PRODUCT CODE	DESCRIPTION	UNITS PER CASE
536070	7.0 mm Replacement Tube with Obturator / Dilator and Inner Cannula	1
536080	8.0 mm Replacement Tube with Obturator / Dilator and Inner Cannula	1
536090	9.0 mm Replacement Tube with Obturator / Dilator and Inner Cannula	1



## Inner Cannula for Per-fit™ Percutaneous Tracheostomy Tubes

PRODUCT CODE	DESCRIPTION	INNER CANNULA I.D.	COLOR CODE	UNITS PER CASE
537070	Inner Cannula for 7.0 mm Tube	6.0 mm	Green	20
537080	Inner Cannula for 8.0 mm Tube	7.0 mm	White	20
537090	Inner Cannula for 9.0 mm Tube	8.0 mm	Blue	20





# Tracheostomy

## Portex® Tracheostomy Tubes (PVC)



<b>Adult Tracheostomy Tubes</b>	<b>p7</b>
<b>Disposable Inner Cannula (D.I.C.®) Tubes</b>	<b>p7</b>
Flex D.I.C.® Tracheostomy Tubes	p7
Cuffed Flex D.I.C.® Tracheostomy Tubes	p7
Uncuffed Flex D.I.C.® Tracheostomy Tubes	p7
D.I.C.® Tracheostomy Tubes	p7
Cuffed D.I.C.® Tracheostomy Tubes	p7
Uncuffed D.I.C.® Tracheostomy Tubes	p7
Cuffed Fenestrated D.I.C.® Tracheostomy Tubes	p7
Uncuffed Fenestrated D.I.C.® Tracheostomy Tubes	p8
Inner Cannula for D.I.C.® and Flex D.I.C.® Tracheostomy Tubes	p8
<b>Laryngectomy Tubes</b>	<b>p8</b>
D.I.C.® Laryngectomy Tubes	p8
Inner Cannula for D.I.C.® Laryngectomy Tubes	p8
<b>Portex® Blue Line ULTRA® Tracheostomy Tubes</b>	<b>p9</b>
Portex® Blue Line ULTRA® Cuffed Tracheostomy Tubes	p9
Portex® Blue Line ULTRA® Cuffed Change Kits with Tube	p9
Portex® Blue Line ULTRA® Uncuffed Tracheostomy Tubes	p9
Portex® Blue Line ULTRA® Uncuffed Change Kits with Tube	p9
Portex® Blue Line ULTRA® Cuffed Fenestrated Tracheostomy Tubes	p9
Portex® Blue Line ULTRA® Uncuffed Fenestrated Tracheostomy Tubes	p10
Portex® Blue Line ULTRA® Suctionaid® Tracheostomy Tubes	p10
Portex® Blue Line ULTRA® Suctionaid® Change Kits with Tube	p10
Inner Cannula for Portex® Blue Line ULTRA® Tracheostomy Tubes	p10

Continued on next page



## Tracheostomy

### Portex® Tracheostomy Tubes (PVC)

<b>Portex® Blue Line® Tracheostomy Tubes</b>	<b>p11</b>
Blue Line® Cuffed Tracheostomy Tubes	p11
Blue Line® Uncuffed Tracheostomy Tubes	p11
Blue Line® Cuffed Fenestrated Tracheostomy Tubes	p11
<b>Talk Tube</b>	<b>p11</b>
Blue Line® Cuffed Trach-Talk™ Tracheostomy Tubes	p11
Inner Cannula for Blue Line® and Blue Line® Trach-Talk™ Tracheostomy Tubes	p11
<b>Adjustable and Extra Length Tracheostomy Tubes</b>	<b>p12</b>
Blue Line® Extra Horizontal Length Adjustable Neck Flange Cuffed Tracheostomy Tubes	p12
Blue Line® Extra Horizontal Length Adjustable Neck Flange Uncuffed Tracheostomy Tubes	p12
Blue Line® Extra Vertical Length Tracheostomy Tubes with Double Cuff	p12
<b>Neonatal and Pediatric Tracheostomy Tubes</b>	<b>p13</b>
Portex® Neonatal Tracheostomy Tubes with 15 mm Connector	p13
Portex® Pediatric Tracheostomy Tubes with 15 mm Connector	p13
Portex® Pediatric Tracheostomy Tubes without Connector	p13





505080

### Cuffed Flex D.I.C.® Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	COLOR CODE	UNITS PER CASE
505060	6.0 mm	6.0 mm	8.2 mm	64.0 mm	Orange	1
505070	7.0 mm	7.0 mm	9.6 mm	70.0 mm	Green	1
505080	8.0 mm	8.0 mm	10.9 mm	73.0 mm	White	1
505090	9.0 mm	9.0 mm	12.3 mm	79.0 mm	Blue	1
505100	10.0 mm	10.0 mm	13.7 mm	79.0 mm	Yellow	1



504080

### Uncuffed Flex D.I.C.® Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	COLOR CODE	UNITS PER CASE
504060	6.0 mm	6.0 mm	8.2 mm	64.0 mm	Orange	1
504070	7.0 mm	7.0 mm	9.6 mm	70.0 mm	Green	1
504080	8.0 mm	8.0 mm	10.9 mm	73.0 mm	White	1
504090	9.0 mm	9.0 mm	12.3 mm	79.0 mm	Blue	1
504100	10.0 mm	10.0 mm	13.7 mm	79.0 mm	Yellow	1



503080

### Cuffed D.I.C.® Tracheostomy Tubes


PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	COLOR CODE	UNITS PER CASE
503060	6.0 mm	6.0 mm	8.5 mm	64.0 mm	Orange	1
503070	7.0 mm	7.0 mm	9.9 mm	70.0 mm	Green	1
503080	8.0 mm	8.0 mm	11.3 mm	73.0 mm	White	1
503090	9.0 mm	9.0 mm	12.6 mm	79.0 mm	Blue	1
503100	10.0 mm	10.0 mm	14.0 mm	79.0 mm	Yellow	1



502080

### Uncuffed D.I.C.® Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	COLOR CODE	UNITS PER CASE
502060	6.0 mm	6.0 mm	8.5 mm	64.0 mm	Orange	1
502070	7.0 mm	7.0 mm	9.9 mm	70.0 mm	Green	1
502080	8.0 mm	8.0 mm	11.3 mm	73.0 mm	White	1
502090	9.0 mm	9.0 mm	12.6 mm	79.0 mm	Blue	1
502100	10.0 mm	10.0 mm	14.0 mm	79.0 mm	Yellow	1



513080

### Cuffed Fenestrated D.I.C.® Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	COLOR CODE	UNITS PER CASE
513060	6.0 mm	6.0 mm	8.5 mm	64.0 mm	Orange	1
513070	7.0 mm	7.0 mm	9.9 mm	70.0 mm	Green	1
513080	8.0 mm	8.0 mm	11.3 mm	73.0 mm	White	1
513090	9.0 mm	9.0 mm	12.6 mm	79.0 mm	Blue	1
513100	10.0 mm	10.0 mm	14.0 mm	79.0 mm	Yellow	1



512080

### Uncuffed Fenestrated D.I.C.<sup>®</sup> Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	COLOR CODE	UNITS PER CASE
512060	6.0 mm	6.0 mm	8.5 mm	64.0 mm	Orange	1
512070	7.0 mm	7.0 mm	9.9 mm	70.0 mm	Green	1
512080	8.0 mm	8.0 mm	11.3 mm	73.0 mm	White	1
512090	9.0 mm	9.0 mm	12.6 mm	79.0 mm	Blue	1
512100	10.0 mm	10.0 mm	14.0 mm	79.0 mm	Yellow	1



### Inner Cannula for D.I.C.<sup>®</sup> and Flex D.I.C.<sup>®</sup> Tracheostomy Tubes

PRODUCT CODE	DESCRIPTION	INNER CANNULA I.D.	COLOR CODE	UNITS PER CASE
526060	Inner Cannula for 6.0 mm Tube	5.0 mm	Orange	60
526070	Inner Cannula for 7.0 mm Tube	6.0 mm	Green	60
526080	Inner Cannula for 8.0 mm Tube	7.0 mm	White	60
526090	Inner Cannula for 9.0 mm Tube	8.0 mm	Blue	60
526100	Inner Cannula for 10.0 mm Tube	9.0 mm	Yellow	60



562080

### D.I.C.<sup>®</sup> Laryngectomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	COLOR CODE	UNITS PER CASE
562080	8.0 mm	8.0 mm	11.3 mm	50.0 mm	White	1
562090	9.0 mm	9.0 mm	12.6 mm	50.0 mm	Blue	1
562100	10.0 mm	10.0 mm	14.0 mm	50.0 mm	Yellow	1



### Inner Cannula for D.I.C.<sup>®</sup> Laryngectomy Tubes

PRODUCT CODE	DESCRIPTION	INNER CANNULA I.D.	COLOR CODE	UNITS PER CASE
563080	Inner Cannula for 8.0 mm Tube	7.0 mm	White	10
563090	Inner Cannula for 9.0 mm Tube	8.0 mm	Blue	10
563100	Inner Cannula for 10.0 mm Tube	9.0 mm	Yellow	10



100/815/080

### Blue Line ULTRA® Cuffed Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
100/815/060	6.0 mm	6.0 mm	9.2 mm	64.5 mm	1
100/815/070	7.0 mm	7.0 mm	10.5 mm	70.0 mm	1
100/815/080	8.0 mm	8.0 mm	11.9 mm	75.5 mm	1
100/815/090	9.0 mm	9.0 mm	13.3 mm	81.0 mm	1
100/815/100	10.0 mm	10.0 mm	14.0 mm	87.5 mm	1

Includes inner cannula with stop ring, disconnection wedge, tracheostomy tube holder, cleaning brush, and patient label.



100/825/080

### Blue Line ULTRA® Cuffed Change Kits with Tube

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
100/825/060	6.0 mm	6.0 mm	9.2 mm	64.5 mm	1
100/825/070	7.0 mm	7.0 mm	10.5 mm	70.0 mm	1
100/825/080	8.0 mm	8.0 mm	11.9 mm	75.5 mm	1
100/825/090	9.0 mm	9.0 mm	13.3 mm	81.0 mm	1
100/825/100	10.0 mm	10.0 mm	14.0 mm	87.5 mm	1

Includes inner cannula with stop ring, guidewire, disconnection wedge, tracheostomy tube holder, cleaning brush, and patient label.



100/816/080

### Blue Line ULTRA® Uncuffed Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
100/816/060	6.0 mm	6.0 mm	9.2 mm	64.5 mm	1
100/816/070	7.0 mm	7.0 mm	10.5 mm	70.0 mm	1
100/816/080	8.0 mm	8.0 mm	11.9 mm	75.5 mm	1
100/816/090	9.0 mm	9.0 mm	13.3 mm	81.0 mm	1
100/816/100	10.0 mm	10.0 mm	14.0 mm	87.5 mm	1

Includes inner cannula with stop ring, disconnection wedge, tracheostomy tube holder, cleaning brush, and patient label.



100/826/080

### Blue Line ULTRA® Uncuffed Change Kits with Tube

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
100/826/060	6.0 mm	6.0 mm	9.2 mm	64.5 mm	1
100/826/070	7.0 mm	7.0 mm	10.5 mm	70.0 mm	1
100/826/080	8.0 mm	8.0 mm	11.9 mm	75.5 mm	1
100/826/090	9.0 mm	9.0 mm	13.3 mm	81.0 mm	1
100/826/100	10.0 mm	10.0 mm	14.0 mm	87.5 mm	1

Includes inner cannula with stop ring, guidewire, disconnection wedge, tracheostomy tube holder, cleaning brush, and patient label.



100/817/080

### Blue Line ULTRA® Cuffed Fenestrated Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
100/817/060	6.0 mm	6.0 mm	9.2 mm	64.5 mm	1
100/817/070	7.0 mm	7.0 mm	10.5 mm	70.0 mm	1
100/817/080	8.0 mm	8.0 mm	11.9 mm	75.5 mm	1
100/817/090	9.0 mm	9.0 mm	13.3 mm	81.0 mm	1
100/817/100	10.0 mm	10.0 mm	14.0 mm	87.5 mm	1

Includes inner cannula with stop ring, disconnection wedge, tracheostomy tube holder, cleaning brush, and patient label.



100/818/080

### Blue Line ULTRA® Uncuffed Fenestrated Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
100/818/060	6.0 mm	6.0 mm	9.2 mm	64.5 mm	1
100/818/070	7.0 mm	7.0 mm	10.5 mm	70.0 mm	1
100/818/080	8.0 mm	8.0 mm	11.9 mm	75.5 mm	1
100/818/090	9.0 mm	9.0 mm	13.3 mm	81.0 mm	1
100/818/100	10.0 mm	10.0 mm	14.0 mm	87.5 mm	1

Includes inner cannula with stop ring, disconnection wedge, tracheostomy tube holder, cleaning brush, and patient label.



100/875/080

### Blue Line ULTRA® Suctionaid® Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
100/875/060	6.0 mm	6.0 mm	9.2 mm	64.5 mm	1
100/875/070	7.0 mm	7.0 mm	10.5 mm	70.0 mm	1
100/875/080	8.0 mm	8.0 mm	11.9 mm	75.5 mm	1
100/875/090	9.0 mm	9.0 mm	13.3 mm	81.0 mm	1
100/875/100	10.0 mm	10.0 mm	14.0 mm	87.5 mm	1

Includes inner cannula with stop ring, disconnection wedge, tracheostomy tube holder, cleaning brush, and patient label.



100/885/080

### Blue Line ULTRA® Suctionaid® Change Kits with Tube

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
100/885/060	6.0 mm	6.0 mm	9.2 mm	64.5 mm	1
100/885/070	7.0 mm	7.0 mm	10.5 mm	70.0 mm	1
100/885/080	8.0 mm	8.0 mm	11.9 mm	75.5 mm	1
100/885/090	9.0 mm	9.0 mm	13.3 mm	81.0 mm	1
100/885/100	10.0 mm	10.0 mm	14.0 mm	87.5 mm	1

Includes inner cannula with stop ring, guidewire, disconnection wedge, tracheostomy tube holder, cleaning brush, and patient label.



100/856/080

### Inner Cannula for Blue Line ULTRA® Tracheostomy Tubes

PRODUCT CODE	DESCRIPTION	INNER CANNULA I.D.	UNITS PER CASE
100/856/060	Inner Cannula for 6.0 mm Tube	5.0 mm	20
100/856/070	Inner Cannula for 7.0 mm Tube	6.0 mm	20
100/856/080	Inner Cannula for 8.0 mm Tube	7.0 mm	20
100/856/090	Inner Cannula for 9.0 mm Tube	8.0 mm	20
100/856/100	Inner cannula for 10.0 mm tube	9.0 mm	20

100/858/060	Inner Cannula for 6.0 mm Tube	5.0 mm	50
100/858/070	Inner Cannula for 7.0 mm Tube	6.0 mm	50
100/858/080	Inner Cannula for 8.0 mm Tube	7.0 mm	50
100/858/090	Inner Cannula for 9.0 mm Tube	8.0 mm	50
100/858/100	Inner Cannula for 10.0 mm Tube	9.0 mm	50



100/782/080

### Blue Line® Cuffed Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
100/782/060	6.0 mm	6.0 mm	8.3 mm	55.0 mm	1
100/782/070	7.0 mm	7.0 mm	9.7 mm	65.0 mm	1
100/782/080	8.0 mm	8.0 mm	11.0 mm	76.0 mm	1
100/782/090	9.0 mm	9.0 mm	12.4 mm	87.0 mm	1
100/782/100	10.0 mm	10.0 mm	13.8 mm	98.0 mm	1



100/781/080

### Blue Line® Uncuffed Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
100/781/060	6.0 mm	6.0 mm	8.3 mm	55.0 mm	1
100/781/070	7.0 mm	7.0 mm	9.7 mm	65.0 mm	1
100/781/080	8.0 mm	8.0 mm	11.0 mm	76.0 mm	1
100/781/090	9.0 mm	9.0 mm	12.4 mm	87.0 mm	1
100/781/100	10.0 mm	10.0 mm	13.8 mm	98.0 mm	1



100/783/080

### Blue Line® Cuffed Fenestrated Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
100/783/070	7.0 mm	7.0 mm	9.7 mm	65.0 mm	1
100/783/080	8.0 mm	8.0 mm	11.0 mm	76.0 mm	1
100/783/090	9.0 mm	9.0 mm	12.4 mm	87.0 mm	1



100/785/080

### Blue Line® Cuffed Trach-Talk™ Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
100/785/070	7.0 mm	7.0 mm	9.7 mm	65.0 mm	1
100/785/080	8.0 mm	8.0 mm	11.0 mm	76.0 mm	1
100/785/090	9.0 mm	9.0 mm	12.4 mm	87.0 mm	1



100/790/080

### Inner Cannula for Blue Line® and Blue Line® Trach-Talk™ Tracheostomy Tubes

PRODUCT CODE	DESCRIPTION	INNER CANNULA I.D.	UNITS PER CASE
100/790/070	Inner Cannula for 7.0 mm Tube	5.0 mm	10
100/790/080	Inner Cannula for 8.0 mm Tube	6.0 mm	10
100/790/090	Inner Cannula for 9.0 mm Tube	7.0 mm	10



100/527/080

## Blue Line® Extra Horizontal Length Adjustable Neck Flange Cuffed Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	VERTICAL LENGTH	ADJUSTABLE HORIZONTAL LENGTH	UNITS PER CASE
100/527/060	6.0 mm	6.0 mm	8.3 mm	23.0 mm	11-31 mm	2
100/527/070	7.0 mm	7.0 mm	9.6 mm	26.0 mm	16-36 mm	2
100/527/080	8.0 mm	8.0 mm	11.0 mm	35.0 mm	17-37 mm	2
100/527/090	9.0 mm	9.0 mm	12.4 mm	38.0 mm	22-42 mm	2
100/527/100	10.0 mm	10.0 mm	13.8 mm	41.0 mm	26-46 mm	2



100/529/080

## Blue Line® Extra Horizontal Length Adjustable Neck Flange Uncuffed Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	VERTICAL LENGTH	ADJUSTABLE HORIZONTAL LENGTH	UNITS PER CASE
100/529/060	6.0 mm	6.0 mm	8.3 mm	23.0 mm	11-31 mm	2
100/529/070	7.0 mm	7.0 mm	9.6 mm	26.0 mm	16-36 mm	2
100/529/080	8.0 mm	8.0 mm	11.0 mm	35.0 mm	17-37 mm	2
100/529/090	9.0 mm	9.0 mm	12.4 mm	38.0 mm	22-42 mm	2
100/529/100	10.0 mm	10.0 mm	13.8 mm	41.0 mm	26-46 mm	2



100/784/080

## Blue Line® Extra Vertical Length Tracheostomy Tubes with Double Cuff

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	VERTICAL LENGTH	TOTAL LENGTH	UNITS PER CASE
100/784/070	7.0 mm	7.0 mm	9.7 mm	41.0 mm	82.0 mm	1
100/784/080	8.0 mm	8.0 mm	11.0 mm	45.0 mm	93.0 mm	1
100/784/090	9.0 mm	9.0 mm	12.4 mm	48.0 mm	102.0 mm	1
100/784/100	10.0 mm	10.0 mm	13.7 mm	52.0 mm	113.0 mm	1



553030

### Portex® Neonatal Tracheostomy Tubes with 15 mm Connector

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
553025	2.5 mm	2.5 mm	4.5 mm	30.0 mm	1
553030	3.0 mm	3.0 mm	5.2 mm	32.0 mm	1
553035	3.5 mm	3.5 mm	5.8 mm	34.0 mm	1



555050

### Portex® Pediatric Tracheostomy Tubes with 15 mm Connector

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
555025	2.5 mm	2.5 mm	4.5 mm	30.0 mm	1
555030	3.0 mm	3.0 mm	5.2 mm	36.0 mm	1
555035	3.5 mm	3.5 mm	5.8 mm	40.0 mm	1
555040	4.0 mm	4.0 mm	6.6 mm	44.0 mm	1
555045	4.5 mm	4.5 mm	7.1 mm	48.0 mm	1
555050	5.0 mm	5.0 mm	7.7 mm	50.0 mm	1
555055	5.5 mm	5.5 mm	8.3 mm	52.0 mm	1



559040

### Portex® Pediatric Tracheostomy Tubes without Connector

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
559025	2.5 mm	2.5 mm	4.5 mm	30.0 mm	1
559030	3.0 mm	3.0 mm	5.2 mm	36.0 mm	1
559035	3.5 mm	3.5 mm	5.8 mm	40.0 mm	1
559040	4.0 mm	4.0 mm	6.6 mm	44.0 mm	1
559045	4.5 mm	4.5 mm	7.1 mm	48.0 mm	1
559050	5.0 mm	5.0 mm	7.7 mm	50.0 mm	1
559055	5.5 mm	5.5 mm	8.3 mm	52.0 mm	1





## Tracheostomy

### Bivona® Tracheostomy Tubes (Silicone)



<b>Adult Silicone Tracheostomy Tubes</b>	<b>p17</b>
<b>Adjustable and Extra Length Tracheostomy Tubes</b>	<b>p17</b>
Bivona® Mid-Range Aire-Cuf® Adjustable Neck Flange Hyperflex™ Tracheostomy Tubes	p17
Bivona® TTS™ Adjustable Neck Flange Hyperflex™ Tracheostomy Tubes	p17
Bivona® Cuffed Fixed Neck Flange Hyperflex™ Tracheostomy Tubes	p17
Bivona® Uncuffed Fixed Neck Flange Hyperflex™ Tracheostomy Tubes	p17
<b>Silicone Tracheostomy Tubes</b>	<b>p17</b>
Bivona® Uncuffed Tracheostomy Tubes	p17
Bivona® TTS™ Cuffed Tracheostomy Tubes	p18
Bivona® TTS™ Fixed Neck Flange Hyperflex™ Tracheostomy Tubes	p18
Bivona® Mid-Range Aire-Cuf® Tracheostomy Tubes	p18
Bivona® Fome-Cuf® Tracheostomy Tubes	p18
Bivona® Sleep Apnea Tracheostomy Tubes	p18
Bivona® Laryngectomy Tubes	p19
<b>Talk Tubes</b>	<b>p19</b>
Bivona® Mid-Range Aire-Cuf® Tracheostomy Tubes with Talk Attachment	p19
Bivona® Fome-Cuf® Tracheostomy Tubes with Talk Attachment	p19
<b>Neonatal and Pediatric Silicone Tracheostomy Tubes</b>	<b>p19</b>
<b>Uncuffed Neonatal and Pediatric Tracheostomy Tubes</b>	<b>p19</b>
Bivona® Uncuffed Neonatal Tracheostomy Tubes	p19
Bivona® Uncuffed Pediatric Tracheostomy Tubes	p19
Bivona® Uncuffed Neonatal FlexTend™ Plus Tracheostomy Tubes	p20
Bivona® Uncuffed Pediatric FlexTend™ Plus Standard Tracheostomy Tubes	p20
Bivona® Uncuffed Pediatric FlexTend™ Plus Extra Length Tracheostomy Tubes	p20

Continued on next page

## Tracheostomy Bivona (Silicone)



<b>Cuffed Neonatal and Pediatric Tracheostomy Tubes</b>	<b>p20</b>
Bivona® TTS™ Cuffed Neonatal Tracheostomy Tubes	p20
Bivona® TTS™ Cuffed Pediatric Tracheostomy Tubes	p20
Bivona® Aire-Cuf® Neonatal Tracheostomy Tubes	p21
Bivona® Aire-Cuf® Pediatric Tracheostomy Tubes	p21
Bivona® Fome-Cuf® Neonatal Tracheostomy Tubes	p21
Bivona® Fome-Cuf® Pediatric Tracheostomy Tubes	p21
<b>Specialty Pediatric Tracheostomy Tubes</b>	<b>p21</b>
Bivona® Uncuffed Adjustable Neck Flange Hyperflex™ Pediatric Tracheostomy Tubes	p21
<b>Customized Tracheostomy Tube Services</b>	<b>p22</b>



75HA80

### Bivona® Mid-Range Aire-Cuf® Adjustable Neck Flange Hyperflex™ Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
75HA60	6.0 mm	6.0 mm	8.7 mm	110.0 mm	1
75HA70	7.0 mm	7.0 mm	10.0 mm	120.0 mm	1
75HA80	8.0 mm	8.0 mm	11.0 mm	130.0 mm	1
75HA90	9.0 mm	9.0 mm	12.3 mm	140.0 mm	1



67HA80

### Bivona® TTS™ Adjustable Neck Flange Hyperflex™ Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
67HA60	6.0 mm	6.0 mm	8.7 mm	110.0 mm	1
67HA70	7.0 mm	7.0 mm	10.0 mm	120.0 mm	1
67HA80	8.0 mm	8.0 mm	11.0 mm	130.0 mm	1
67HA90	9.0 mm	9.0 mm	12.3 mm	140.0 mm	1



75FHXL70

### Bivona® Cuffed Extra Length Fixed Flange Hyperflex™ Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
75FHXL60	6.0 mm	6.0 mm	8.7 mm	100.0 mm	1
75FHXL70	7.0 mm	7.0 mm	10.0 mm	110.0 mm	1
75FHXL80	8.0 mm	8.0 mm	11.0 mm	120.0 mm	1
75FHXL90	9.0 mm	9.0 mm	12.3 mm	130.0 mm	1



60AFHXL80

### Bivona® Uncuffed Extra Length Fixed Flange Hyperflex™ Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
60AFHXL60	6.0 mm	6.0 mm	8.7 mm	100.0 mm	1
60AFHXL70	7.0 mm	7.0 mm	10.0 mm	110.0 mm	1
60AFHXL80	8.0 mm	8.0 mm	11.0 mm	120.0 mm	1
60AFHXL90	9.0 mm	9.0 mm	12.3 mm	130.0 mm	1



60A180

### Bivona® Uncuffed Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
60A150	5.0 mm	5.0 mm	7.3 mm	60.0 mm	1
60A160	6.0 mm	6.0 mm	8.7 mm	70.0 mm	1
60A170	7.0 mm	7.0 mm	10.0 mm	80.0 mm	1
60A180	8.0 mm	8.0 mm	11.0 mm	88.0 mm	1
60A190	9.0 mm	9.0 mm	12.3 mm	98.0 mm	1
60A195	9.5 mm	9.5 mm	13.3 mm	98.0 mm	1



670180

### Bivona® TTS™ Cuffed Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
670150	5.0 mm	5.0 mm	7.3 mm	60.0 mm	1
670160	6.0 mm	6.0 mm	8.7 mm	70.0 mm	1
670165	6.5 mm	6.5 mm	9.4 mm	70.0 mm	1
670170	7.0 mm	7.0 mm	10.0 mm	80.0 mm	1
670175	7.5 mm	7.5 mm	10.4 mm	80.0 mm	1
670180	8.0 mm	8.0 mm	11.0 mm	88.0 mm	1
670185	8.5 mm	8.5 mm	11.8 mm	88.0 mm	1
670190	9.0 mm	9.0 mm	12.3 mm	98.0 mm	1
670195	9.5 mm	9.5 mm	13.3 mm	98.0 mm	1



67FH70

### Bivona® TTS™ Fixed Neck Flange Hyperflex™ Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
67FH60	6.0 mm	6.0 mm	8.7 mm	70.0 mm	1
67FH70	7.0 mm	7.0 mm	10.0 mm	80.0 mm	1
67FH80	8.0 mm	8.0 mm	11.0 mm	88.0 mm	1
67FH90	9.0 mm	9.0 mm	12.3 mm	98.0 mm	1



750180

### Bivona® Mid-Range Aire-Cuf® Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
750150	5.0 mm	5.0 mm	7.3 mm	60.0 mm	1
750160	6.0 mm	6.0 mm	8.7 mm	70.0 mm	1
750170	7.0 mm	7.0 mm	10.0 mm	80.0 mm	1
750180	8.0 mm	8.0 mm	11.0 mm	88.0 mm	1
750190	9.0 mm	9.0 mm	12.3 mm	98.0 mm	1
750195	9.5 mm	9.5 mm	13.3 mm	98.0 mm	1



850180

### Bivona® Fome-Cuf® Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
850150	5.0 mm	5.0 mm	7.3 mm	60.0 mm	1
850160	6.0 mm	6.0 mm	8.7 mm	70.0 mm	1
850170	7.0 mm	7.0 mm	10.0 mm	80.0 mm	1
850180	8.0 mm	8.0 mm	11.0 mm	88.0 mm	1
850190	9.0 mm	9.0 mm	12.3 mm	98.0 mm	1
850195	9.5 mm	9.5 mm	13.3 mm	98.0 mm	1



SAT656

### Bivona® Sleep Apnea Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
SAT605	6.0 mm	6.0 mm	9.0 mm	68.0 mm	1
SAT656	6.5 mm	6.5 mm	10.0 mm	73.0 mm	1
SAT757	7.5 mm	7.5 mm	11.0 mm	79.0 mm	1



BOSL5L

### Bivona® Laryngectomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
BOSL1L	9.5 mm	9.5 mm	12.0 mm	85.0 mm	1
BOSL1S	9.5 mm	9.5 mm	12.0 mm	55.0 mm	1
BOSL3L	11.0 mm	11.0 mm	14.0 mm	105.0 mm	1
BOSL3S	11.0 mm	11.0 mm	14.0 mm	55.0 mm	1
BOSL5L	13.0 mm	13.0 mm	16.0 mm	85.0 mm	1
BOSL5S	13.0 mm	13.0 mm	16.0 mm	55.0 mm	1



755180

### Bivona® Mid-Range Aire-Cuf® Tracheostomy Tubes with Talk Attachment

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
755150	5.0 mm	5.0 mm	7.3 mm	60.0 mm	1
755160	6.0 mm	6.0 mm	8.7 mm	70.0 mm	1
755170	7.0 mm	7.0 mm	10.0 mm	80.0 mm	1
755180	8.0 mm	8.0 mm	11.0 mm	88.0 mm	1
755190	9.0 mm	9.0 mm	12.3 mm	98.0 mm	1
755195	9.5 mm	9.5 mm	13.3 mm	98.0 mm	1



855180

### Bivona® Fome-Cuf® Tracheostomy Tubes with Talk Attachment

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
855150	5.0 mm	5.0 mm	7.3 mm	60.0 mm	1
855160	6.0 mm	6.0 mm	8.7 mm	70.0 mm	1
855170	7.0 mm	7.0 mm	10.0 mm	80.0 mm	1
855180	8.0 mm	8.0 mm	11.0 mm	88.0 mm	1
855190	9.0 mm	9.0 mm	12.3 mm	98.0 mm	1
855195	9.5 mm	9.5 mm	13.3 mm	98.0 mm	1



60N025

### Bivona® Uncuffed Neonatal Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
60N025	2.5 mm	2.5 mm	4.0 mm	30.0 mm	1
60N030	3.0 mm	3.0 mm	4.7 mm	32.0 mm	1
60N035	3.5 mm	3.5 mm	5.3 mm	34.0 mm	1
60N040	4.0 mm	4.0 mm	6.0 mm	36.0 mm	1



60P045

### Bivona® Uncuffed Pediatric Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
60P025	2.5 mm	2.5 mm	4.0 mm	38.0 mm	1
60P030	3.0 mm	3.0 mm	4.7 mm	39.0 mm	1
60P035	3.5 mm	3.5 mm	5.3 mm	40.0 mm	1
60P040	4.0 mm	4.0 mm	6.0 mm	41.0 mm	1
60P045	4.5 mm	4.5 mm	6.7 mm	42.0 mm	1
60P050	5.0 mm	5.0 mm	7.3 mm	44.0 mm	1
60P055	5.5 mm	5.5 mm	8.0 mm	46.0 mm	1



60NFP30

### Bivona® Uncuffed Neonatal FlexTend™ Plus Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	PROXIMAL LENGTH	DISTAL LENGTH	UNITS PER CASE
60NFP25	2.5 mm	2.5 mm	4.0 mm	15.0 mm	30.0 mm	1
60NFP30	3.0 mm	3.0 mm	4.7 mm	17.0 mm	32.0 mm	1
60NFP35	3.5 mm	3.5 mm	5.3 mm	20.0 mm	34.0 mm	1
60NFP40	4.0 mm	4.0 mm	6.0 mm	22.0 mm	36.0 mm	1



60PFS45

### Bivona® Uncuffed Pediatric FlexTend™ Plus Standard Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	PROXIMAL LENGTH	DISTAL LENGTH	UNITS PER CASE
60PFS25	2.5 mm	2.5 mm	4.0 mm	20.0 mm	38.0 mm	1
60PFS30	3.0 mm	3.0 mm	4.7 mm	20.0 mm	39.0 mm	1
60PFS35	3.5 mm	3.5 mm	5.3 mm	25.0 mm	40.0 mm	1
60PFS40	4.0 mm	4.0 mm	6.0 mm	25.0 mm	41.0 mm	1
60PFS45	4.5 mm	4.5 mm	6.7 mm	30.0 mm	42.0 mm	1
60PFS50	5.0 mm	5.0 mm	7.3 mm	30.0 mm	44.0 mm	1
60PFS55	5.5 mm	5.5 mm	8.0 mm	30.0 mm	46.0 mm	1



60PFP50

### Bivona® Uncuffed Pediatric FlexTend™ Plus Extra Length Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	PROXIMAL LENGTH	DISTAL LENGTH	UNITS PER CASE
60PFP35	3.5 mm	3.5 mm	5.3 mm	38.0 mm	40.0 mm	1
60PFP40	4.0 mm	4.0 mm	6.0 mm	40.0 mm	44.0 mm	1
60PFP45	4.5 mm	4.5 mm	6.7 mm	43.0 mm	48.0 mm	1
60PFP50	5.0 mm	5.0 mm	7.3 mm	45.0 mm	50.0 mm	1
60PFP55	5.5 mm	5.5 mm	8.0 mm	48.0 mm	52.0 mm	1



67N025

### Bivona® TTS™ Cuffed Neonatal Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
67N025	2.5 mm	2.5 mm	4.0 mm	30.0 mm	1
67N030	3.0 mm	3.0 mm	4.7 mm	32.0 mm	1
67N035	3.5 mm	3.5 mm	5.3 mm	34.0 mm	1
67N040	4.0 mm	4.0 mm	6.0 mm	36.0 mm	1



67P040

### Bivona® TTS™ Cuffed Pediatric Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
67P025	2.5 mm	2.5 mm	4.0 mm	38.0 mm	1
67P030	3.0 mm	3.0 mm	4.7 mm	39.0 mm	1
67P035	3.5 mm	3.5 mm	5.3 mm	40.0 mm	1
67P040	4.0 mm	4.0 mm	6.0 mm	41.0 mm	1
67P045	4.5 mm	4.5 mm	6.7 mm	42.0 mm	1
67P050	5.0 mm	5.0 mm	7.3 mm	44.0 mm	1
67P055	5.5 mm	5.5 mm	8.0 mm	46.0 mm	1



65N025

### Bivona® Aire-Cuf® Neonatal Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
65N025	2.5 mm	2.5 mm	4.0 mm	30.0 mm	1
65N030	3.0 mm	3.0 mm	4.7 mm	32.0 mm	1
65N035	3.5 mm	3.5 mm	5.3 mm	34.0 mm	1
65N040	4.0 mm	4.0 mm	6.0 mm	36.0 mm	1



65P050

### Bivona® Aire-Cuf® Pediatric Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
65P025	2.5 mm	2.5 mm	4.0 mm	38.0 mm	1
65P030	3.0 mm	3.0 mm	4.7 mm	39.0 mm	1
65P035	3.5 mm	3.5 mm	5.3 mm	40.0 mm	1
65P040	4.0 mm	4.0 mm	6.0 mm	41.0 mm	1
65P045	4.5 mm	4.5 mm	6.7 mm	42.0 mm	1
65P050	5.0 mm	5.0 mm	7.3 mm	44.0 mm	1
65P055	5.5 mm	5.5 mm	8.0 mm	46.0 mm	1



85N030

### Bivona® Fome-Cuf® Neonatal Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
85N025	2.5 mm	2.5 mm	4.0 mm	30.0 mm	1
85N030	3.0 mm	3.0 mm	4.7 mm	32.0 mm	1
85N035	3.5 mm	3.5 mm	5.3 mm	34.0 mm	1
85N040	4.0 mm	4.0 mm	6.0 mm	36.0 mm	1



85P035

### Bivona® Fome-Cuf® Pediatric Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	TUBE LENGTH	UNITS PER CASE
85P025	2.5 mm	2.5 mm	4.0 mm	38.0 mm	1
85P030	3.0 mm	3.0 mm	4.7 mm	39.0 mm	1
85P035	3.5 mm	3.5 mm	5.3 mm	40.0 mm	1
85P040	4.0 mm	4.0 mm	6.0 mm	41.0 mm	1
85P045	4.5 mm	4.5 mm	6.7 mm	42.0 mm	1
85P050	5.0 mm	5.0 mm	7.3 mm	44.0 mm	1
85P055	5.5 mm	5.5 mm	8.0 mm	46.0 mm	1



60HA35

### Bivona® Uncuffed Adjustable Neck Flange Hyperflex™ Pediatric Tracheostomy Tubes

PRODUCT CODE	TUBE SIZE	TUBE I.D.	TUBE O.D.	MAXIMUM USABLE TUBE LENGTH	UNITS PER CASE
60HA25	2.5 mm	2.5 mm	4.0 mm	55.0 mm	1
60HA30	3.0 mm	3.0 mm	4.7 mm	60.0 mm	1
60HA35	3.5 mm	3.5 mm	5.3 mm	65.0 mm	1
60HA40	4.0 mm	4.0 mm	6.0 mm	70.0 mm	1
60HA45	4.5 mm	4.5 mm	6.7 mm	75.0 mm	1
60HA50	5.0 mm	5.0 mm	7.3 mm	80.0 mm	1
60HA55	5.5 mm	5.5 mm	8.0 mm	85.0 mm	1



## Tracheostomy Bivona® Customized Tubes

### Customized Tracheostomy Tube Services

Special problems require special solutions. At Smiths Medical we understand the necessity of putting your patient's needs first. We offer a comprehensive range of components and component styles from our broad silicone tube line to accommodate your special tracheostomy tube requirements.

Together we review your patient's requirements, provide you with convenient order templates, and upon your written approval, we will assemble a tracheostomy tube to address your patient needs.

Our "turn-around" goal from receipt of your signed template to shipment is seven days for non-sterile tubes that can be steam autoclave. To obtain priority service and assure the fastest possible delivery, ask your Customer Service Representative about our "Same-Day" and "Bivona Express" Services.

Our broad range of customization options:

#### Tube Shaft Style

- Standard Silicone
- Hyperflex™ Wire Reinforced Silicone
- Fenestrations
- Stoma Cushions- Stomainer shafts

#### Curvature

- Standard Anatomical Pre-curve Shaft
- Hyperflex™ Wire Reinforced Flexible
- Anatomically Adapting Shaft

#### Length

- Standard Extensive Range

#### I.D. / O.D.

- Comprehensive range to choose from

#### Cuff Design

- Fome-Cuf®
- Aire-Cuf®
- TTS™
- Cuffless
- Multiple Cuff Configurations
- Stoma Cushion Cuff

#### Cuff Position

- To your specification upon request

#### Neck Flange

- Adult
- Neonatal ("V" and Straight)
- Pediatric ("V" and Straight)
- Fixed
- Adjustable










## Tracheostomy Bivona® Tracheostoma Vents

04

### Tracheostoma Vents

PRODUCT CODE	SIZE	I.D.	O.D.	LENGTH	UNITS PER CASE
<b>BV0818</b>	8	9.5 mm	12.0 mm	18.0 mm	1
<b>BV0827</b>	8	9.5 mm	12.0 mm	27.0 mm	1
<b>BV0836</b>	8	9.5 mm	12.0 mm	36.0 mm	1
<b>BV0855</b>	8	9.5 mm	12.0 mm	55.0 mm	1
<b>BV0918</b>	9	10.5 mm	13.5 mm	18.0 mm	1
<b>BV0927</b>	9	10.5 mm	13.5 mm	27.0 mm	1
<b>BV0936</b>	9	10.5 mm	13.5 mm	36.0 mm	1
<b>BV0955</b>	9	10.5 mm	13.5 mm	55.0 mm	1
<b>BV1018</b>	10	12.0 mm	15.0 mm	18.0 mm	1
<b>BV1027</b>	10	12.0 mm	15.0 mm	27.0 mm	1
<b>BV1036</b>	10	12.0 mm	15.0 mm	36.0 mm	1
<b>BV1055</b>	10	12.0 mm	15.0 mm	55.0 mm	1
<b>BV1218</b>	12	13.5 mm	17.0 mm	18.0 mm	1
<b>BV1227</b>	12	13.5 mm	17.0 mm	27.0 mm	1
<b>BV1236</b>	12	13.5 mm	17.0 mm	36.0 mm	1
<b>BV1255</b>	12	13.5 mm	17.0 mm	55.0 mm	1





## Tracheostomy Cricothyrotomy Kits

### Mini-Trach® II Cricothyrotomy Kit

PRODUCT CODE	DESCRIPTION	UNITS PER CASE
500100	Mini-Trach® II Cricothyrotomy Kit	10



### Nu-Trake® Adult Emergency Cricothyrotomy Kit

PRODUCT CODE	DESCRIPTION	UNITS PER CASE
B10100	Nu-Trake® Cricothyrotomy Kit	1



# Tracheostomy

## Tracheostomy Ancillary Products and Accessories



<b>Velcro® Tracheostomy Tube Holders</b>	<b>p30</b>
<b>Tracheostomy Tube Straps</b>	<b>p30</b>
<b>Tracheostomy Care Trays</b>	<b>p30</b>
<b>Tracheostomy Accessories</b>	<b>p30</b>
<b>Tracheostomy Educational Aids</b>	<b>p31</b>



## Velcro® Tracheostomy Tube Holders

PRODUCT CODE	DESCRIPTION	UNITS PER CASE
520000	Adult	100
520001	Pediatric	100
520002	Infant	100



## Tracheostomy Tube Straps

PRODUCT CODE	DESCRIPTION	MAXIMUM NECK SIZE	UNITS PER CASE
69-0600TS	Adult Standard One Piece, 1½" Wide Extra Cushion	22"	100
69-0675TS	Adjustable Child to Adult, 1½" Wide Extra Cushion Two Components / Cut-to-Fit	24"	100
69-0680TS	Adjustable Child to Adult, 1" Wide Extra Cushion Two Components / Cut-to-Fit	18"	100
69-0685TS	Hook and Loop Closure	20"	100
69-0650	Replacement Tracheostomy Cushion Pads		250



6950

## Tracheostomy Care Trays

PRODUCT CODE	DESCRIPTION	UNITS PER CASE
6950	Full Component Tracheostomy Care Tray with Plastic Forceps	20
6970	Full Component Tracheostomy Care Tray without Plastic Forceps	20
6952	Economy Tracheostomy Care Tray	20
6956	6.0 mm D.I.C.® Tracheostomy Care Tray	20
6957	7.0 mm D.I.C.® Tracheostomy Care Tray	20
6958	8.0 mm D.I.C.® Tracheostomy Care Tray	20
6959	9.0 mm D.I.C.® Tracheostomy Care Tray	20
6960	10.0 mm D.I.C.® Tracheostomy Care Tray	20

## Tracheostomy Educational Aids

PRODUCT CODE	DESCRIPTION	UNITS PER CASE
LT-2151	Portex® Per-fit™ Percutaneous Tracheostomy Kit Clinical Guide	1
FCVIDEO	Bivona® Fome-Cuf® Video Training Kit	1
MDTRAK	Model Trachea Kit	1
NTVIDEO	Bivona® Nu-Trake® Video Training Tape	1



## Tracheostomy Accessories

PRODUCT CODE	DESCRIPTION	UNITS PER CASE
519000	Portex® Decannulation Cap	10
REDCAP	Bivona® Decannulation Cap	26
CMD006	Cuff Maintenance Device with 60 cc Syringe and Stopcock	6
CON15S	SidePort™ Connectors	50
DRYPORT	DryPort™ Kit	5
DWG300	15 mm Disconnect Wedge	25
525350	Portex® PEEP-Keep™ Dual Axis Swivel Adapter (non-sterile)	50
525351	Portex® PEEP-Keep™ Dual Axis Swivel Adapter (sterile)	50
CARSWL	Carden Swivel Connector	1
COSTF1	Saf T Flo® Connector with Port	10
COSTF2	Saf T Flo® Connector without Port	10
002400	Portex® First Breath® Adult Tracheostomy Mask	50
002405	Portex® First Breath® Pediatric Tracheostomy Mask	50



519000



REDCAP



CON15S



CMD006



DRYPORT



525350



CARSWL

# Index

REORDER NO.	PAGE	REORDER NO.	PAGE	REORDER NO.	PAGE	REORDER NO.	PAGE
002400	31	100/815/100	9	100/885/090	10	535070	2
002405	31	100/816/060	9	100/885/100	10	535080	2
100/527/060	12	100/816/070	9	500100	27	535090	2
100/527/070	12	100/816/080	9	502060	7	536070	2
100/527/080	12	100/816/090	9	502070	7	536080	2
100/527/090	12	100/816/100	9	502080	7	536090	2
100/527/100	12	100/817/060	9	502090	7	537070	2
100/529/060	12	100/817/070	9	502100	7	537080	2
100/529/070	12	100/817/080	9	503060	7	537090	2
100/529/080	12	100/817/090	9	503070	7	553025	13
100/529/090	12	100/817/100	9	503080	7	553030	13
100/529/100	12	100/818/060	10	503090	7	553035	13
100/572/000	2	100/818/070	10	503100	7	555025	13
100/572/LIT	2	100/818/080	10	504060	7	555030	13
100/591/070	2	100/818/090	10	504070	7	555035	13
100/591/080	2	100/818/100	10	504080	7	555040	13
100/591/090	2	100/825/060	9	504090	7	555045	13
100/781/060	11	100/825/070	9	504100	7	555050	13
100/781/070	11	100/825/080	9	505060	7	555055	13
100/781/080	11	100/825/090	9	505070	7	559025	13
100/781/090	11	100/825/100	9	505080	7	559030	13
100/781/100	11	100/826/060	9	505090	7	559035	13
100/782/060	11	100/826/070	9	505100	7	559040	13
100/782/070	11	100/826/080	9	512060	8	559045	13
100/782/080	11	100/826/090	9	512070	8	559050	13
100/782/090	11	100/826/100	9	512080	8	559055	13
100/782/100	11	100/856/060	10	512090	8	562080	8
100/783/070	11	100/856/070	10	512100	8	562090	8
100/783/080	11	100/856/080	10	513060	7	562100	8
100/783/090	11	100/856/090	10	513070	7	563080	8
100/784/070	12	100/856/100	10	513080	7	563090	8
100/784/080	12	100/858/060	10	513090	7	563100	8
100/784/090	12	100/858/070	10	513100	7	60A150	17
100/784/100	12	100/858/080	10	519000	31	60A160	17
100/785/070	11	100/858/090	10	520000	30	60A170	17
100/785/080	11	100/858/100	10	520001	30	60A180	17
100/785/090	11	100/875/060	10	520002	30	60A190	17
100/790/070	11	100/875/070	10	525350	31	60A195	17
100/790/080	11	100/875/080	10	525351	31	60AFHXL60	17
100/790/090	11	100/875/090	10	526060	8	60AFHXL70	17
100/815/060	9	100/875/100	10	526070	8	60AFHXL80	17
100/815/070	9	100/885/060	10	526080	8	60AFHXL90	17
100/815/080	9	100/885/070	10	526090	8	60HA25	21
100/815/090	9	100/885/080	10	526100	8	60HA30	21

# Index

REORDER NO.	PAGE	REORDER NO.	PAGE	REORDER NO.	PAGE
60HA35	21	670160	18	750190	18
60HA40	21	670165	18	750195	18
60HA45	21	670170	18	755150	19
60HA50	21	670175	18	755160	19
60HA55	21	670180	18	755170	19
60N025	19	670185	18	755180	19
60N030	19	670190	18	755190	19
60N035	19	670195	18	755195	19
60N040	19	67FH60	18	75FHXL60	17
60NFP25	20	67FH70	18	75FHXL70	17
60NFP30	20	67FH80	18	75FHXL80	17
60NFP35	20	67FH90	18	75FHXL90	17
60NFP40	20	67HA60	17	75HA60	17
60P025	19	67HA70	17	75HA70	17
60P030	19	67HA80	17	75HA80	17
60P035	19	67HA90	17	75HA90	17
60P040	19	67N025	20	850150	18
60P045	19	67N030	20	850160	18
60P050	19	67N035	20	850170	18
60P055	19	67N040	20	850180	18
60PFP35	20	67P025	20	850190	18
60PFP40	20	67P030	20	850195	18
60PFP45	20	67P035	20	855150	19
60PFP50	20	67P040	20	855160	19
60PFP55	20	67P045	20	855170	19
60PFS25	20	67P050	20	855180	19
60PFS30	20	67P055	20	855190	19
60PFS35	20	6950	30	855195	19
60PFS40	20	6952	30	85N025	21
60PFS45	20	6956	30	85N030	21
60PFS50	20	6957	30	85N035	21
60PFS55	20	6958	30	85N040	21
65N025	21	6959	30	85P025	21
65N030	21	6960	30	85P030	21
65N035	21	6970	30	85P035	21
65N040	21	69-0600TS	30	85P040	21
65P025	21	69-0650	30	85P045	21
65P030	21	69-0675TS	30	85P050	21
65P035	21	69-0680TS	30	85P055	21
65P040	21	69-0685TS	30	B10100	27
65P045	21	750150	18	BOSL1L	19
65P050	21	750160	18	BOSL1S	19
65P055	21	750170	18	BOSL3L	19
670150	18	750180	18	BOSL3S	19
				BOSL5L	19
				BOSL5S	19
				BV0818	25
				BV0827	25
				BV0836	25
				BV0855	25
				BV0918	25
				BV0927	25
				BV0936	25
				BV0955	25
				BV1018	25
				BV1027	25
				BV1036	25
				BV1055	25
				BV1218	25
				BV1227	25
				BV1236	25
				BV1255	25
				CARSWL	31
				CMD006	31
				CON15S	31
				COSTF1	31
				COSTF2	31
				DRYPORT	31
				DWG300	31
				FCVIDEO	30
				LT-2151	30
				MDTRAK	30
				NTVIDEO	30
				REDCAP	31
				SAT605	18
				SAT656	18
				SAT757	18

THE DETAILS GIVEN IN THIS LEAFLET ARE CORRECT AT THE TIME OF GOING TO PRESS. THE COMPANY RESERVES THE RIGHT TO IMPROVE THE EQUIPMENT SHOWN.

**Smiths Medical ASD, Inc.**

160 Weymouth Street, Rockland, MA 02370

Phone: (800) 848-1757 Fax: (800) 621-2590

[www.smiths-medical.com](http://www.smiths-medical.com)

**Smiths Medical, part of the global technology business Smiths Group**

Portex, Bivona, ULTRApert, Per-fit, D.I.C., Blue Line ULTRA, Suctionaid, Blue Line, Trach-Talk, TTS, Hyperflex, Aire-Cuf, Fome-Cuf, FlexTend, Mini-Trach, Nu-Trake, Side-Port, Dry-Port, PEEP-Keep, Saf T Flo, First Breath and Smiths design mark are trademarks of the Smiths Medical family of companies. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners. Please refer to Instructions For Use supplied with the product for detailed instructions, warnings, and cautions. © 2007 Smiths Medical family of companies. All rights reserved. SM118 8.2007



**smiths medical**