



Introducing the new SureTemp® Plus 690

Now taking temperatures is easier, safer, and more accurate.
And taking care of all your patients is more convenient.

WelchAllyn®

Advancing Frontline Care™



- **Easy to use, intuitive operation**
- **Removable probe well reduces the risk of cross-contamination**
- **Improved accuracy for even more reliable readings**
- **Greater versatility for use with all your patients, all the time**

The New SureTemp® Plus 690—
with the speed and accuracy
you require, and the versatility
you need for use with all your
patients, all the time.

When we set out to improve the Welch Allyn SureTemp, we knew we had our work cut out for us. After all, when you start with the thermometer that helped set the standard for speed and accuracy, the expectations run pretty high.

So how did we do? See for yourself by taking a look at the innovative new SureTemp Plus 690— easier, safer and more accurate than the original, and still built to withstand constant use in demanding environments.

Easy: an easy-to-read LCD display with flashing icons that clearly communicate in any language, and an intuitive user interface and ergonomic shape make using the SureTemp Plus 690 a breeze.

Safe: patient safety is improved thanks to an interchangeable probe well that reduces the risk of cross-contamination.

Fast and Accurate: in addition to improved accuracy, the SureTemp 690 is one of the fastest thermometers on the market, providing 4-second oral, 10-second pediatric axillary, 15-second adult axillary, and 10-second rectal measurements.

Versatile: to top it all off, you can now use the SureTemp Plus 690 with all your patients, from infants to adults— whether you need to take oral, axillary or rectal temperatures.

Interchangeable probe with push-button probe cover ejection

Removable probe well for easy cleaning

Convenient storage for 25 probe covers


Large, easy-to-read LCD

C/F button for toggling between Fahrenheit and Celsius

Prominent recall button for display of last temperature taken

Ergonomic design with rubberized grips

Optional wall mount has storage space for 50 additional probe covers



Battery life indicator

User selectable icon and/or words for identification of thermometer mode

Mode selection button for selecting oral, pediatric axillary, or adult axillary



Interchangeable rectal probe and well reduces risk of cross-contamination

The new SureTemp® Plus 690

WelchAllyn®

Advancing Frontline Care™

The New

SureTemp® Plus 690—

setting the standard of

care even higher



SureTemp® Plus 690 specifications

Calibration Accuracy.....	±0.2°F(±0.1°C) (Monitor Mode)
Temperature Predict Time:	
Oral.....	Approx. 4 secs.
Adult Axillary (18 years and older).....	Approx. 15 secs.
Pediatric Axillary (17 years and younger).....	Approx. 10 secs.
Rectal.....	Approx. 10 secs.
Unit Dimensions.....	8.46" x 3.18" x 2.43" (215mm x 81mm x 62mm)
Weight.....	12.6 ounces (357 grams)
Ambient Temperature Range.....	50.0°F(10°C) - 104°F (40.0°C)
Patient Temperature Range.....	80.0°F(26.7°C) - 110°F (43.3°C)
Power Source.....	(3)1.5Vdc AA batteries
Battery Life.....	Approx. 6000 temperatures
Warranty.....	2 year limited instrument, 1 year limited probe, and 90 day limited probe well

ordering information

thermometers

- 01690-200** SureTemp Plus Electronic Thermometer (Model 690), Oral Probe with Oral Probe Well, 2 Year Limited Warranty
- 01690-201** SureTemp Plus Electronic Thermometer (Model 690), Rectal Probe with Rectal Probe Well, 2 Year Limited Warranty
- 01690-300** SureTemp Plus Electronic Thermometer (Model 690), One Per Room, Wall Mount, 9ft. Oral Probe with Oral Probe Well, 2 Year Limited Warranty (*must be locked to the wall*)
- 01690-301** SureTemp Plus Electronic Thermometer (Model 690), One Per Room, Wall Mount, 9ft. Rectal Probe with Rectal Probe Well, 2 Year Limited Warranty (*must be locked to the wall*)
- 01690-400** SureTemp Plus Electronic Thermometer (Model 690), Wall Mount, Oral Probe with Oral Probe Well, 2 Year Limited Warranty
- 01690-700** SureTemp Plus Electronic Thermometer (Model 690), Wall Mount, 9ft. Oral Probe with Oral Probe Well, 2 Year Limited Warranty, with Rolling Stand and Holder for Spare Probe & Well

accessories

- 21326-0000** Wall Holder
- 21330-0000** Holder for Spare Probe & Well
- 06138-0000** Calibration Key Assy, 690/692
- 02892-0000** Probe & Well Kit 4FT Rectal
- 02892-1000** Probe & Well Kit 9FT Rectal
- 02893-0000** Probe & Well Kit 4FT Oral
- 02893-1000** Probe & Well Kit 9FT Oral
- 02891-0000** Probe Well, Blue
- 02891-1000** Probe Well, Red
- 08273-0000** Thermometer Stand with Holder for Spare Probe & Well
- 21333-0000** Bracket, M600 Stand Upgrade, with Holder for Spare Probe & Well
- 21331-0000** Bracket Assembly, M690/692 to Tycos Stand
- 98024-0000** Extended Warranty M690

replacement probe covers

- 05031-101** Case of 1,000 probe covers

WelchAllyn®

Advancing Frontline Care™

Call Welch Allyn to order the SureTemp® Plus 690
7420 Carroll Road, San Diego, CA 92121-2334 USA
Tel: 800.535.6663 / 858.621.6600 • Fax: 858.621.6611
www.welchallyn.com

SureTemp is a registered trademark of Welch Allyn, Inc. © 2003

SM2546 Rev. B