



## Introducing the new SureTemp® Plus 692

Speed and accuracy you require.  
Safety and security you demand.

**WelchAllyn®**

Advancing Frontline Care™



- **Easy to use, intuitive operation**
- **Removable probe well reduces the risk of cross-contamination**
- **Improved accuracy for even more reliable readings**
- **Innovative security features that cut down on the chance of theft**

## The New SureTemp® Plus 692— with the speed and accuracy you require, and the safety and security you demand

When we set out to improve the Welch Allyn SureTemp, we knew we had our work cut out for us. After all, when you start with the thermometer that helped set the standard for speed and accuracy, the expectations run pretty high.


So how did we do? See for yourself by taking a look at the innovative new SureTemp Plus 692—easier, safer and more accurate than the original, and still built to withstand constant use in demanding environments.

**Easy:** an easy-to-read LCD display with flashing icons that clearly communicate in any language, and an intuitive user interface and ergonomic shape make using the SureTemp Plus 692 a breeze.

**Safe:** patient safety is improved thanks to an interchangeable probe well that reduces the risk of cross-contamination.

**Fast and Accurate:** in addition to improved accuracy, the SureTemp 692 is one of the fastest thermometers on the market, providing 4-second oral, 10-second pediatric axillary, 15-second adult axillary, and 10-second rectal measurements.

**Secure:** to top it all off, our new security features make sure the unit stays where it is supposed to be. A location ID on the display and other innovative features help prevent theft.



Interchangeable probe with push-button probe cover ejection

Removable probe well for easy cleaning

Convenient storage for 25 probe covers

C/F button for toggling between Fahrenheit and Celsius

ID Field for electronic labeling of thermometer



Battery life indicator

User selectable icon and/or words for identification of thermometer mode

Large, easy-to-read, backlit LCD

Mode selection button for selecting oral, pediatric axillary, or adult axillary

Pulse timer button for timing pulse and respiration

Ergonomic design with rubberized grips

Prominent recall button for display of last temperature taken

Storage space for 50 additional probe covers in the wall mount (Wall mount standard with the 692 model)

Optional rolling stand for security and portability

Interchangeable rectal probe and well reduces risk of cross-contamination



# The new SureTemp® Plus 692

**WelchAllyn®**

Advancing Frontline Care™

The New

SureTemp® Plus 692—

setting the standard of  
care even higher



## SureTemp® Plus 692 specifications

Calibration Accuracy.....	±0.2°F(±0.1°C) (Monitor Mode)
Temperature Predict Time:	
Oral.....	Approx. 4 secs.
Adult Axillary (18 years and older).....	Approx. 15 secs.
Pediatric Axillary (17 years and younger).....	Approx. 10 secs.
Rectal.....	Approx. 10 secs.
Unit Dimensions.....	8.46" x 3.18" x 2.43" (215mm x 81mm x 62mm)
Weight.....	12.6 ounces (357 grams)
Ambient Temperature Range.....	50.0°F(10°C) - 104°F (40.0°C)
Patient Temperature Range.....	80.0°F(26.7°C) - 110°F (43.3°C)
Power Source.....	(3)1.5Vdc AA batteries
Battery Life.....	Approx. 5000 temperatures
Warranty.....	3 year limited instrument, 1 year limited probe, and 90 day limited probe well

## ordering information

### thermometers

- 01692-200** SureTemp Plus Electronic Thermometer (Model 692), Wall Mount, Security System with ID Location Field, and Oral Probe with Oral Probe Well, 3 Year Limited Warranty
- 01692-201** SureTemp Plus Electronic Thermometer (Model 692), Wall Mount, Security System with ID Location Field, and Rectal Probe with Rectal Probe Well, 3 Year Limited Warranty
- 01692-300** SureTemp Plus Electronic Thermometer (Model 692), One Per Room, Wall Mount, 9ft. Oral Probe with Oral Probe Well, 3 Year Limited Warranty (*must be locked to the wall*)
- 01692-301** SureTemp Plus Electronic Thermometer (Model 692), One Per Room, Wall Mount, 9ft. Rectal Probe with Rectal Probe Well, 3 Year Limited Warranty (*must be locked to the wall*)
- 01692-700** SureTemp Plus Electronic Thermometer (Model 692), Wall Mount, 9ft. Oral Probe with Oral Probe Well, 3 Year Limited Warranty, with Rolling Stand and Holder for Spare Probe & Well

### accessories

- 21326-0000** Wall Holder
- 21330-0000** Holder for Spare Probe & Well
- 06138-0000** Calibration Key Assy, 690/692
- 02892-0000** Probe & Well Kit 4FT Rectal
- 02892-1000** Probe & Well Kit 9FT Rectal
- 02893-0000** Probe & Well Kit 4FT Oral
- 02893-1000** Probe & Well Kit 9FT Oral
- 02891-0000** Probe Well, Blue
- 02891-1000** Probe Well, Red
- 08273-0000** Thermometer Stand with Holder for Spare Probe & Well
- 21333-0000** Bracket, M600 Stand Upgrade, with Holder for Spare Probe & Well
- 21331-0000** Bracket Assembly, M690/692 to Tyco Stand
- 98025-0000** Extended Warranty M692

### replacement probe covers

- 05031-101** Case of 1,000 probe covers

**WelchAllyn®**

Advancing Frontline Care™

Call Welch Allyn to order the SureTemp® Plus 692  
7420 Carroll Road, San Diego, CA 92121-2334 USA  
Tel: 800.535.6663 / 858.621.6600 • Fax: 858.621.6611  
[www.welchallyn.com](http://www.welchallyn.com)

SureTemp is a registered trademark of Welch Allyn, Inc. © 2003

SM2550 Rev. B