


# Dimensions


Aurelia® Dimensions™ Fitted Powder Free  
Latex Examination Gloves



**Hand-Specific Comfort  
Superb Ergonomics**



- Left and right hand precision-fit gloves for most comfortable fit
- Reduces risk of carpal tunnel syndrome
- Fingertip-textured finish for precise gripping
- Polymer-coated helps gloves slide on with ease, keeping hand smooth



**AURELIA**  
Sign of Comfort™


[www.aureliagloves.com](http://www.aureliagloves.com) | 1-877-AURELIA (1-877-287-3542)

# Dimensions

## Product Description

**Aurelia® Dimensions™ Fitted**

Powder Free Latex Examination Gloves are made from 100% natural Malaysian rubber, Hand-Specific, creamy white color, Fingertip textured and non-sterile. Aurelia gloves conform and are tested to the highest international standards. For more information visit [www.aureliagloves.com](http://www.aureliagloves.com)



## Aurelia® Dimensions™ Specifications

### TENSILE TEST RESULTS (unaged)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	30.0	970
Small	31.1	970
Medium	31.0	970
Large	31.2	960
Extra Large	31.3	970
ASTM 3578:05	18 MPa minimum	650% minimum

### AGED TENSILE TEST RESULTS (7 Days at 70° C)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	28.0	900
Small	28.0	920
Medium	27.5	910
Large	27.3	890
Extra Large	27.9	900
ASTM 3578:05	14 MPa minimum	500% minimum

## DIMENSIONS

	Length (inches/mm)	Thickness (Single-wall)	
		Palm (mils/mm)	Finger (mils/mm)
<b>6.5</b> Aurelia Dimensions ASTM 3578:05	9.4 / 240 8.7 / 220 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
<b>7</b> Aurelia Dimensions ASTM 3578:05	9.4 / 240 8.7 / 220 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
<b>7.5</b> Aurelia Dimensions ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
<b>8</b> Aurelia Dimensions ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
<b>8.5</b> Aurelia Dimensions ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum

### ITEM NUMBER

6.5	7	7.5	8	8.5
69925	69926	69927	69928	69929

## Ordering Information

To order from your preferred distributor, specify size and item number. Packing: 100 pieces per dispenser box, 10 dispenser boxes per shipping carton.



Authorized Aurelia Distributor:

**CAUTION STATEMENT:** This product contains natural rubber latex which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.

**AURELIA Div. Supermax Inc.**

2225 White Oak Circle, Aurora, IL 60502 | Toll Free: 1-877-AURELIA | Tel: 1-630-898-8886 | Fax: 1-630-898-8855  
Web: [www.aureliagloves.com](http://www.aureliagloves.com) | Email: [info@aureliagloves.com](mailto:info@aureliagloves.com)